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EN-Nesse nasal aspirator

General safety rules

- This product is intended for use **only under adult supervision** and in accordance with these instructions.
- Before first and each subsequent use, **wash all parts that come into contact with the nose and mouth** according to the instructions.
- Before each use, **check the condition of tips and adapters** - if they are damaged, cracked or deformed, **do not use them**.
- This product is not a toy - keep out of reach of children when not in use.
- Keep the packaging in case you need to consult.

Using an aspirator for infants and children

- Use the aspirator only to **remove nasal secretions** gently and **briefly** - each suction should not last long so as not to irritate the mucous membrane.
- Do not insert the tip too deep into your nose; apply it only to the entrance of your nostrils.

- If your child shows discomfort, nosebleeds or has signs of inflammation, **discontinue use** and contact your doctor.
- To facilitate removal of thick secretions, consider using **saline solution** before aspiration (as per pediatric practices).

Hygiene and maintenance

- After each use , **disassemble and wash the components** that come into contact with secretions, preferably in warm water with a mild detergent; dry before reassembling and storing.
- The model is compatible with Gledlig , Hyggelig and Apring breast pumps , **clean the adapters** and all connecting elements after each use to prevent the transfer of bacteria.
- Store all parts in **a clean and dry place** , avoiding moisture and direct sunlight.

Warnings and restrictions on use

- The product **is not intended for use other than removing nasal secretions** .
- Do not use the aspirator in **your ears, mouth or eyes** .
- Do not use if your child has **wounds, inflammation, infections or nosebleeds** without consulting a doctor.
- If using the oral aspirator function, **avoid direct vacuum suction stronger than the recommended level** to avoid causing discomfort or damage to the nasal/oral mucosa.
- If you notice **any unexpected vibrations, sounds or unusual signals from the product** , discontinue use.

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